

Model: Dipping tank for work benches sheet steel, 97 l

Order number: 4.273.22

Correct use: Safety containers and cleaning containers



capacity: approx. 97 litres outside dimensions: 895 x 496 x 336 mm (W x D x H) basin (inside): 889 x 406 x 336 mm (W x D x H) For positioning on a workbench. Basin, lid and frame made from galvanised sheet steel, red coated. self-closing lid in case of fire "3/4"- drain plug in basin bottom for selective emptying Testing: FM-approved

Category: Cleaning tanks Outer height [mm]: 305, Outer width [mm]: 895, Outer depth [mm]: 406, Gross weight [kg]: 25,00, Net weight [kg]: 22,00,