

# **Model: Steel dipping tank capacity 19 l**

## **Order number: 4.276.05**

Correct use: Safety containers and cleaning containers



Dim: approx. 360 x 360 (dia x H in mm), outside: powder coated, red, inside: tin-plated, Self-closing lid in case of fire via safety fuse. Markings prepared according to GefStoffV and BGV A8 testing: FM/ UL- tested

Category: Cleaning containers Outer height [mm]: 360, Gross weight [kg]: 5,00, Net weight [kg]: 4,10,