

Model: Sheet steel dipping tank with PE-insert, capacity 19 l Order number: 4.276.15

Correct use: Safety containers and cleaning containers



Dim: 360 x 360 (diaxH in mm) outside: powder coated, red, inside: tin-plated, Self-closing lid in case of fire via safety fuse. Markings prepared according to GefStoffV and BGV A8 testing: FM/ UL-tested

Category: Cleaning containers Outer height [mm]: 360, Gross weight [kg]: 5,50, Net weight [kg]: 5,00,