

Model: Workbench collection container with foot pedal, capacity 20 l

Order number: 5.901.00

Correct use: Safety containers and cleaning containers



Dim: 290 x 397 (dia x H mm), sheet steel, exterior red powder coated, interior tin-plated, self-closing lid using foot pedal for opening, Markings prepared according to GefStoffV and BGV A8 testing: FM/UL-approved

Category: for solid hazardous substances Outer height [mm]: 397, Gross weight [kg]: 4,50, Net weight [kg]: 4,10, Outer diameter [mm]: 397,