

## **Model: Workbench collection container with hand opening**

### **Order number: 5.901.10**

Correct use: Safety containers and cleaning containers



Capacity 20 l, Dim: 290 x 397 (dia x H mm), sheet steel, exterior red powder coated, interior tin-plated, self-closing lid, opening by hand, Markings prepared according to GefStoffV and BGV A8 testing: GS-sign and FM/UL-tested

Category: for solid hazardous substances Outer height [mm]: 397, Gross weight [kg]: 4,00, Net weight [kg]: 3,70,