

Model: Workbench collection container table-top version, capacity 8 l Order number: 5.902.00

Correct use: Safety containers and cleaning containers



Dim.: 245 x 230 (dia x H in mm), sheet steel, exterior red powder coated, interior tin-plated, without perforated floor edge, with side handles, self-closing lid, opening by hand, markings prepared according to GefStoffV and BGV A8. testing: FM-tested

Category: for solid hazardous substances Outer height [mm]: 230, Gross weight [kg]: 3,00, Net weight [kg]: 3,00,