

Model: Workbench collection container with foot pedal, capacity 50 l

Order number: 5.905.00

Correct use: Safety containers and cleaning containers



Dim: 393 x 505 (dia x H mm), sheet steel, exterior red powder coated, interior tin-plated, self-closing lid, using foot pedal for opening, Markings prepared according to GefStoffV and BGV A8 testing: FM/UL-tested

Category: for solid hazardous substances Outer height [mm]: 505, Gross weight [kg]: 7,50, Net weight [kg]: 6,80,