

Model: Workbench collection container with foot pedal, capacity 80 l
Order number: 5.907.00

Correct use: Safety containers and cleaning containers



Dim:approx. 458 x 597 (dia x H mm), sheet steel, exterior red powder coated, interior tin-plated, self-closing lid, using foot pedal for opening, Markings prepared according to GefStoffV and BGV A8 testing: FM/UL-tested

Category:	for solid hazardous substances
Outer height [mm]:	597
Net weight [kg]:	10,00
Outer diameter [mm]:	458