

Model: Workbench collection container with foot pedal, capacity 80 l

Order number: 5.907.00

Correct use: Safety containers and cleaning containers



Dim: approx. 458 x 597 (dia x H mm), sheet steel, exterior red powder coated, interior tin-plated, self-closing lid, using foot pedal for opening, Markings prepared according to GefStoffV and BGV A8 testing: FM/UL-tested

Category: for solid hazardous substances Outer height [mm]: 597, Gross weight [kg]: 10,00, Net weight [kg]: 10,00, Outer diameter [mm]: 458,