

**Model: Compressed air media pump explosion protected (set)**  
**Order number: 50-70000-052**

Correct use: Supply (active storage)



"MEMBRANE PUMP 1/4"" max. 6bar air pressure Capacity: 0-16 L/min weight: 3,2kg Parts tangent to liquids made from: Teflon, nickel-plated brass, PP, PA and stainless steel.

Category:	Equipment
Outer width [mm]:	145
Outer depth [mm]:	195
Outer height [mm]:	137
Net weight [kg]:	3,20