

Model: Compressed air media pump explosion protected (set)

Order number: 50-70000-052

Correct use: Supply (active storage)



"MEMBRANE PUMP 1/4"" max. 6bar air pressure Capacity: 0-16 L/min weight: 3,2kg Parts tangent to liquids made from: Teflon, nickel-plated brass, PP, PA and stainless steel.

Category: Equipment Outer height [mm]: 137, Outer width [mm]: 145, Outer depth [mm]: 195, Gross weight [kg]: 4,00, Net weight [kg]: 3,20,