

# Model: Gas cylinder trolley for the transport of a 50 l cylinder

## Order number: 7.00.600

Correct use: Compressed gas cylinders storage



Blue painted steel tube construction. Including securing chain. With solid rubber tyres. Dimensions (W x D x H) approx.: 480 x 580 x 1050 mm

Category: Handling of gas cylinders Outer height [mm]: 580, Outer width [mm]: 1.050, Outer depth [mm]: 480, Gross weight [kg]: 16,00, Net weight [kg]: 11,50,