

# **Model: Gas cylinder trolley, load: 100 kg, for 2 x 10 l cylinders**

## **Order number: 7.00.617**

Correct use: Compressed gas cylinders storage



Suitable for the transport of 2 x 10 litre gas cylinders. With solid rubber tyres. Load bearing capacity: 100 kg. Including toolbox and wrench SW 30/32. Colour: RAL 7015, coarse structure. Size appropriate for the storage in safety cabinets for gas cylinders: models SUPREME plus M,L,XL,XXL as well as models ECO plus M,L,XL,XXL and gas cylinder depots M,L,XL, Dimensions (W x D x H) approx.: 356 x 317 x 1000 mm

Category: Handling of gas cylinders Outer height [mm]: 1.000, Outer width [mm]: 360, Outer depth [mm]: 320, Gross weight [kg]: 11,50, Net weight [kg]: 10,50,