

Model: Gas cylinder trolley, load: 50 kg, for 1 x 10 l cylinder
Order number: 7.00.618

Correct use: Compressed gas cylinders storage



Suitable for the transport of 1 x 10 litre gas cylinder. With solid rubber tyres. Load bearing capacity: 50 kg. Including toolbox and wrench SW 30/32. Colour: RAL 7015, coarse structure. Size appropriate for the storage in safety cabinets for gas cylinders: models SUPREME plus M,L,XL,XXL as well as models ECO plus M,L,XL,XXL and gas cylinder depots M,L,XL, Dimensions (W x D x H) approx.: 275 x 317 x 1000 mm

Category:	Handling of gas cylinders
Outer width [mm]:	280
Outer depth [mm]:	320
Outer height [mm]:	1.000
Net weight [kg]:	10,00