

Model: Gas cylinder trolley for 2 x 40/50 litre cylinders
Order number: 7.00.620

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Suitable for the transport of 2 x 50 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Load bearing capacity: 200 kg. Dimensions (W x D x H) approx.: 770 x 780 x 1350 mm

Category:	Handling of gas cylinders
Outer width [mm]:	1.350
Outer depth [mm]:	770
Outer height [mm]:	780
Net weight [kg]:	31,00