

Model: Gas cylinder trolley for 2 x 40/50 litre cylinders

Order number: 7.00.620

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Suitable for the transport of 2 x 50 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Load bearing capacity: 200 kg. Dimensions (W x D x H) approx.: 770 x 780 x 1350 mm

Category: Handling of gas cylinders Outer height [mm]: 780, Outer width [mm]: 1.350, Outer depth [mm]: 770, Gross weight [kg]: 32,00, Net weight [kg]: 31,00,