

Model: Gas cylinder pallets for 8 x gas cylinder dia. 250 mm Order number: 7.00.810

Correct use: Compressed gas cylinders storage



Surface hot-dip galvanised according to EN ISO1461 dim in mm: 1200 x 860 x 1000 (W x D x H)
Max. load bearing capacity: 700 kg

Category: Handling of gas cylinders Outer height [mm]: 1.000, Outer width [mm]: 1.200, Outer depth [mm]: 860, Gross weight [kg]: 94,00, Net weight [kg]: 94,00,