

Model: Gas cylinder trolley for propane gas cylinders
Order number: 7.01.000

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Maximum number of propane gas cylinders: 1 x 11 kg. Blue painted steel tube construction. Including securing chain. Finger guards for safe grip. Dimensions (W x D x H) approx.: 450 x 550 x 1030 mm

Category:	Handling of gas cylinders
Outer width [mm]:	1.030
Outer depth [mm]:	450
Outer height [mm]:	550
Net weight [kg]:	11,00