

Model: Gas cylinder trolley for propane gas cylinders
Order number: 7.01.010

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Maximum number of propane gas cylinders: 1 x 33 kg. Blue painted steel tube construction. Including securing chain. Red PVC handle for safe grip. Dimensions (W x D x H) approx.: 485 x 550 x 1050 mm

Category:	Handling of gas cylinders
Outer width [mm]:	1.050
Outer depth [mm]:	485
Outer height [mm]:	550
Net weight [kg]:	13,00