

Model: Gas cylinder trolley for propane gas cylinders

Order number: 7.01.010

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Maximum number of propane gas cylinders: 1 x 33 kg. Blue painted steel tube construction. Including securing chain. Red PVC handle for safe grip. Dimensions (W x D x H) approx.: 485 x 550 x 1050 mm

Category: Handling of gas cylinders Outer height [mm]: 550, Outer width [mm]: 1.050, Outer depth [mm]: 485, Gross weight [kg]: 14,00, Net weight [kg]: 13,00,