

Model: Gas cylinder trolley for 2 x 20 litre gas cylinders

Order number: 7.01.020

Correct use: Compressed gas cylinders storage



With solid rubber tyres. Suitable for the transport of 2 x 20 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Red PVC handle for safe grip. Dimensions (W x D x H) approx.: 530 x 500 x 1050 mm

Category: Handling of gas cylinders Outer height [mm]: 500, Outer width [mm]: 1.050, Outer depth [mm]: 530, Gross weight [kg]: 14,00, Net weight [kg]: 13,50,