

Model: Gas cylinder trolley for 2 x 40/50 litre cylinders
Order number: 7.01.030

Correct use: Compressed gas cylinders storage



With solid rubber tyres and fold-out tail wheel. Suitable for the transport of 2 x 50 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Finger guards for safe grip. Dimensions (W x D x H) approx.: 770 x 800 x 1350 mm

Category:	Handling of gas cylinders
Outer width [mm]:	1.350
Outer depth [mm]:	770
Outer height [mm]:	800
Net weight [kg]:	42,00