

# **Model: Gas cylinder trolley for 2 x 40/50 litre cylinders**

## **Order number: 7.01.030**

Correct use: Compressed gas cylinders storage



With solid rubber tyres and fold-out tail wheel. Suitable for the transport of 2 x 50 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Finger guards for safe grip. Dimensions (W x D x H) approx.: 770 x 800 x 1350 mm

Category: Handling of gas cylinders Outer height [mm]: 800, Outer width [mm]: 1.350, Outer depth [mm]: 770, Gross weight [kg]: 43,00, Net weight [kg]: 42,00,