

Model: Gas cylinder trolley for 2 x 40/50 litre cylinders

Order number: 7.01.035

Correct use: Compressed gas cylinders storage



With solid rubber tyres and fold-out tail wheel. Suitable for the transport of 2 x 50 litre gas cylinders. Blue painted steel tube construction. Including securing chain. Finger guards for safe grip. With lockable toolbox, large hose and torch holder and two welding rod holders. Dimensions (W x D x H) approx.: 770 x 800 x 1600 mm

Category: Handling of gas cylinders Outer height [mm]: 800, Outer width [mm]: 1.600, Outer depth [mm]: 770, Gross weight [kg]: 46,00, Net weight [kg]: 45,00,