

**Model: Gas cylinder container with roof**  
**Order number: 7.80.128**

Correct use: Compressed gas cylinders outdoor storage



for 60 gas cylinders dia = 220 mm for outdoor storage of gas cylinders in accordance with the regulations. surface hot-dip galvanised according to EN ISO1461 delivery in modular system with assembly instructions for self-assembly incl. double wing door, W x H 2160 x 1985 mm with profile cylinder dimensions (outside) in mm: 2400 x 1500 x 2160 (WxDxH) Design according to TRGS 510.

Category:	Gas cylinder containers
Outer width [mm]:	2.400
Outer depth [mm]:	1.500
Outer height [mm]:	2.160
Net weight [kg]:	255,00