

Model: Gas cylinder container with roof

Order number: 7.80.128

Correct use: Compressed gas cylinders outdoor storage



for 60 gas cylinders dia = 220 mm for outdoor storage of gas cylinders in accordance with the regulations. surface hot-dip galvanised according to EN ISO1461 delivery in modular system with assembly instructions for self-assembly incl. double wing door, W x H 2160 x 1985 mm with profile cylinder dimensions (outside) in mm: 2400 x 1500 x 2160 (WxDxH) Design according to TRGS 510.

Category: Gas cylinder containers Outer height [mm]: 2.160, Outer width [mm]: 2.400, Outer depth [mm]: 1.500, Gross weight [kg]: 255,00, Net weight [kg]: 255,00,