

Model: Gas cylinder container with roof
Order number: 7.80.154

Correct use: Compressed gas cylinders outdoor storage



for 104 gas cylinders dia = 220 mm for outdoor storage of gas cylinders in accordance with the regulations. surface hot-dip galvanised according to EN ISO1461 delivery in modular system with assembly instructions for self-assembly dimensions (outside) in mm: 3100 x 2100 x 2180 (WxDxH) Design according to TRGS 510.

Category:	Gas cylinder containers
Outer width [mm]:	3.100
Outer depth [mm]:	2.100
Outer height [mm]:	2.180
Net weight [kg]:	320,00