

**Model: Gas cylinder box with 1 door perforated door and side panels**  
**Order number: 7.90.300**

Correct use: Compressed gas cylinders outdoor storage



perforated doors and side panels for natural ventilation zinc coated sheet steel for max. 8 x 33 kg gas cylinders Dim in mm: 840 x 690 x 1500 (W x D x H) with lock.

Category:	Gas cylinder box
Outer width [mm]:	840
Outer depth [mm]:	690
Outer height [mm]:	1.500
Net weight [kg]:	50,00