

Model: Gas cylinder box with 1 door perforated door and side panels

Order number: 7.90.300

Correct use: Compressed gas cylinders outdoor storage



perforated doors and side panels for natural ventilation zinc-coated sheet steel for max. 8 x 33 kg gas cylinders Dim in mm: 840 x 690 x 1500 (W x D x H) with lock.

Category: Gas cylinder box Outer height [mm]: 1.500, Outer width [mm]: 840, Outer depth [mm]: 690, Gross weight [kg]: 82,00, Net weight [kg]: 50,00,