

Model: Gas cylinder box with 2 doors perforated doors and side panels **Order number: 7.90.350**

Correct use: Compressed gas cylinders outdoor storage



perforated doors and side panels for natural ventilation zinc-coated sheet steel for max. 16 x 33 kg gas cylinders Dim in mm: 1680 x 690 x 1500 (W x D x H) with lock.

Category:	Gas cylinder box
Outer width [mm]:	1.680
Outer depth [mm]:	690
Outer height [mm]:	1.500
Net weight [kg]:	98,00