

Model: Gas cylinder box with 2 doors perforated doors and side panels

Order number: 7.90.350

Correct use: Compressed gas cylinders outdoor storage



perforated doors and side panels for natural ventilation zinc coated sheet steel for max. 16 x 33 kg gas cylinders Dim in mm: 1680 x 690 x 1500 (W x D x H) with lock.

Category: Gas cylinder box Outer height [mm]: 1.500, Outer width [mm]: 1.680, Outer depth [mm]: 690, Gross weight [kg]: 129,00, Net weight [kg]: 98,00,