

Model: GASFLASCHEN-DEPOT
Order number: 70-201040-003

Correct use: Compressed gas cylinders outdoor storage



Classical door technology (door handle). For the storage of pressurised gas cylinders outdoors.
Construction and colour:

- Carcass and doors sheet steel, RAL 9002 Fittings and features:
- Cylinder lock for safe access
- Cylinder holder incl 3 retaining belts
- Fitting holder Further technical information: www.dueperthal.com

Category:	Gas cylinder depot - outside storage
Model size	L
Outer width [mm]:	954
Outer depth [mm]:	400
Outer height [mm]:	2.053
Net weight [kg]:	103,00
Number of storage levels	1

Number of storage areas per storage level	1
Door operation	classical technology
Door technology	wing door technology with inspection window