

**Model: GASFLASCHEN-DEPOT**  
**Order number: 70-201340-003**

Correct use: Compressed gas cylinders outdoor storage



Classical door technology (door handle). For the storage of pressurised gas cylinders outdoors.  
Construction and colour:

- Carcass and doors sheet steel, RAL 9002 Fittings and features:
- Cylinder lock for safe access
- Cylinder holder incl 4 retaining belts
- Fitting holder Further technical information: [www.dueperthal.com](http://www.dueperthal.com)

Category:	Gas cylinder depot - outside storage
Model size	XL
Outer width [mm]:	1.254
Outer depth [mm]:	400
Outer height [mm]:	2.053
Net weight [kg]:	125,00
Number of storage levels	1
Number of storage areas per storage level	1

Door operation	classical technology
Door technology	wing door technology with inspection window