

# Model: GASFLASCHEN-DEPOT

## Order number: 70-201340-003

Correct use: Compressed gas cylinders outdoor storage



Classical door technology (door handle). For the storage of pressurised gas cylinders outdoors.  
Construction and colour:

- Carcass and doors sheet steel, RAL 9002 Fittings and features:
- Cylinder lock for safe access
- Cylinder holder incl 4 retaining belts
- Fitting holder Further technical information: [www.dueperthal.com](http://www.dueperthal.com)

Category: Gas cylinder depot - outside storage Model size: XL, Outer height [mm]: 2.053, Outer width [mm]: 1.254, Outer depth [mm]: 400, Gross weight [kg]: 125,00, Net weight [kg]: 125,00, Number of storage levels: 1, Number of storage areas per storage level: 1, Door operation: classical technology, Door technology: wing door technology with inspection window,