

**Model: Elegant hygiene column made of stainless steel**  
**Order number: 9.19.730**

Correct use: Disinfection



Suitable for 1 liter container of disinfectant (9.19.100) Hygiene column for flexible installation with large base for stable standing Removable drip tray as drip protection in fixed mounting Guide frame on the rear side for easy insertion and exchange of disinfectant 1 liter container (arcticle 9.19.100)

Category:	Disinfection
Outer width [mm]:	420
Outer depth [mm]:	290
Outer height [mm]:	1.200
Net weight [kg]:	12,90