

Model: Set of disinfectant for wall mounting with dispenser holder made of stainless steel

Order number: 9.19.820

Correct use: Disinfection



incl. 1 liter container of disinfectant. Set of disinfectant consisting of: Stainless steel dispenser holder with drip protection for wall mounting, mounting material included disinfectant 1 liter container (article 9.19.100) metering pump (article 9.19.710)

Category: Disinfection Outer height [mm]: 205, Outer width [mm]: 94, Outer depth [mm]: 112, Gross weight [kg]: 0,00,