

Model: Set of disinfectant as table stand with dispenser holder

Order number: 9.19.825

Correct use: Disinfection



Made of stainless steel incl. 1 liter container of disinfectant Set of desinfectant consisting of:
Stainless steel dispenser holder with drip protection for wall mounting, desinfectant 1 liter
container (article 9.19.100) metering pump (article 9.19.710)

Category: Disinfection Outer height [mm]: 209, Outer width [mm]: 149, Outer depth [mm]: 149,
Gross weight [kg]: 0,00,