

# Model: Eye rinsing station with two bottles of 500 ml

**Order number: 9.31.805**

Correct use: Emergency showers



sterile sodium chloride solution (0.9%). With wall mounting bracket and integrated pictogram. The station is additionally equipped with a mirror and allows for a particularly correct and quick rinsing at workplaces where a longer rinsing time is required. Dimensions wall station (H x W x D): 290 x 230 x 80 mm. Shelf life: 3 years. Rinsing time: approx. 10 minutes. According to DIN EN 15154-4:2009.

Category: Eye showers Gross weight [kg]: 2,30, Net weight [kg]: 2,00,