

# Model: Eye emergency station with wall mounting bracket

Order number: 9.31.845

Correct use: Emergency showers



With two bottles: 1 x 200 ml of pH-neutral phosphate buffer solution (4.9%) and 1 x 500 ml of sterile sodium chloride solution (0.9%). With integrated pictogram and mirror. Incl. eye emergency plan. Especially suited for workplaces where acids, lye and foreign bodies can get into the eye. Dimensions wall station (H x W x D): 290 x 230 x 80 mm. Shelf life: 3 years. Rinsing time: approx. 7 minutes. According to DIN EN 15154-4:2009.

|                  |             |
|------------------|-------------|
| Category:        | Eye showers |
| Net weight [kg]: | 1,30        |