

Model: Eye rinsing station with two bottles, 500 ml

Order number: 9.32.005

Correct use: Emergency showers



in dust-proof wall box. With two bottles of 500 ml sterile sodium chloride solution (0.9%), integrated pictogram and mirror. Especially suited for workplaces, where dust and dirt accumulate (e.g. workshops and garages). In a first-aid-case very easy and quick to open. Dimensions wall box (H x W x D): 265 x 227 x 110 mm. Shelf life: 3 years. Rinsing time: approx. 10 minutes. According to DIN EN 15154-4:2009.

Category: Eye showers Gross weight [kg]: 2,30, Net weight [kg]: 2,00,