

Model: Eye emergency station in wall box
Order number: 9.32.045

Correct use: Emergency showers



With two bottles: 1 x 200 ml pH-neutral phosphate buffer solution (4.9%) and 1 x 500 ml sterile sodium chloride solution (0.9%). With integrated pictogram and mirror. Incl. eye emergency plan. Especially suited for workplaces where acids, lye and foreign bodies can get into the eye. Dimensions wall station (H x W x D): 265 x 227 x 110 mm. Shelf life: 3 years. Rinsing time: approx. 7 minutes. According to DIN EN 15154-4:2009.

| | |
|------------------|-------------|
| Category: | Eye showers |
| Net weight [kg]: | 1,30 |